Development of robust iPSC-based \alpha-Synuclein, Tau and TDP-43 aggregation models for drug discovery



Jessica Koepke, Isabel Onofre, Carola van Berkel, Tatiana Wust, Leo Smit, Joyson Dcosta, Chiara Formica, Kimberley Riegman, Benoit Samson-Couterie and Shushant Jain Ncardia Services BV Leiden, The Netherlands, support@ncardia.com, www.ncardia.com

Disclosures: None

Background

Proteinopathies are diseases caused by protein misfolding and self-aggregation, which leads to altered neuronal function and neurodegeneration. Examples of proteinopathies include Parkinson's Disease (PD), Alzheimer's Disease (AD) and Amyotrophic lateral sclerosis (ALS).

Proteinopathies are complex diseases difficult to model in vitro and in vivo. Many therapies that showed promise in animal studies have failed in human clinical trials, emphasizing a translation gap in the drug discovery process.

Developing human physiologically relevant disease models is of high importance to identify and validate drugs therapeutic potential with higher confidence of clinical success. Human induced pluripotent stem cells (hiPSCs) have the potential to be differentiated into any cell type, retain patient-specific genetic backgrounds, mimic clinically-relevant human (patho-)physiology and respond appropriately to candidate therapeutics. These characteristics make hiPSCs an excellent tool for drug discovery.

In this study, Ncardia developed three in vitro assays based on hiPSC-derived neurons to model the aggregation of α-synuclein, Tau and TDP-43 – key hallmarks for diseases like Parkinson's, Alzheimer's or Amyotrophic Lateral Sclerosis.

Figure A shows the experimental set-up, timeline and endpoints of Ncardia's AD, PD and ALS pathology model. Assays in purple and in green can be combined and multiplexed





Postulated model for protein aggregation mechanisms in neurodegenerative diseases

- Misfolding of specific characteristic disease related proteins is suggested to be linked to disease progression, resulting in aggregation and fibril formation of these proteins.
- One disease protein aggregates and loses its function or additionally shows a toxic gain of function.

Bahareh Eftekharzadeh, Bradley T. Hyman, Susanne Wegmann Structural studies on the mechanism of protein aggregation in age related neurodegenerative diseases, Mechanisms of Ageing and Development, Volume 156, 2016, Pages 1-13, ISSN 0047-6374, https://doi.org/10.1016/j.mad.2016.03.001.



3. Aggregation of Tau and pTau - Quantifying aggregation

125₇ DAPI-pTau **DAPI-Tau** Merge DAPI-MC-1 DAPI-pTau **DAPI-Tau** Merge DAPI-MC-1 100-



Figure A Tau overexpression by transduction in combination with chronic treatment with Tau PFFs (condition transduced +PFF) showed a notable increase in misfolded Tau (MC-1 in purple) and pTau (green) signals, resembling Tau pathology observed in Alzheimer's disease (magnification 20x)



Figure B Fixation of neurons with methanol, removes the phosphorylated Tau (green) and MC-1 positive Tau (purple) but in condition transduced +PFFs, Tau signal (orange) is still present, indicating the formation of insoluble tau aggregates (magnification 20x).



Figure C Quantification of the normalized parameters (intensity and puncta counts of MC-1 and pTau) of untransduced cells (UNT in purple) compared to cells transduced and treated with PFFs (TR+ PFFs in turquoise), expressed as % to condition UNT set to 0%.

Plotted values are averages ±SD. Statistical significance calculated with unpaired t-test * $p \le 0.05$; ** $p \le 0.01$; *** p≤0.001; **** p≤0000.1

Overexpression of TAU combined with chronic exposure to PFFs induced an Alzheimer's disease phenotype

4. Mis-localization and aggregation of TDP-43 and downregulation STMN2



motor neurons treated with the stressor (magnification 40x).

100%.

Motor neurons chronically treated with a stressor show mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant

* iCell® Motor Neurons 01279 from FUJIFILM Cellular Dynamics, Inc.

consistently and robustly increases TDP-43 aggregation.

Conclusions

condition set to 100%.

• Neuronal co-cultures were used to quantify disease relevant phenotypes for α-synuclein or TAU aggregation as well as the formation intermediate phosphorylated species after treatment with α-synuclein and Tau recombinant preformed fibrils (PFFs) using a high content imaging platform. Inhibitors and activators of protein degradation produced the expected changes in protein aggregation in a concentrationdependent manner. Stressor-treated mutant and wild type iPSC-derived motor neurons (hiPSC-MN) showed disease-specific mis-localization of TDP-43 to the cytoplasm, aggregation of TDP-43, reduction of STMN2 and appearance of the truncated variant.

Figure B Quantification of the normalized parameter (ratio of cytoplasmic/nuclear TDP43) of wild type and mutated

Figure C Quantification of the normalised parameter (TDP-43 aggregation) by HTRF of wild type and mutated TDP-43

(M337V mutation) motor neurons, compared to wild type neurons (in purple) and expressed as % to condition set to 100%. The mutant motor neurons show a suggestive increase in TDP-43 aggregation but the addition of a stressor,

TDP-43 (M337V mutation) motor neurons, compared to wild type neurons (in purple) and expressed as % to

- We have established a suite of robust, clinically relevant in vitro assays (Z-factor >0.5) for the aggregation of α-synuclein, Tau and TDP-43 using human iPSC-derived neuronal subtypes. These assays are performed in a scalable wellplate format and are fully automated to support drug developers at any stage of their discovery process.
- We successfully modelled and evaluated disease-linked phenotypes relevant to AD, PD and ALS, among other neurodegenerative disorders, using complimentary assays. Altogether, offering the opportunity to gain a holistic understanding of the efficacy of therapeutics targeting aggregation.



Questions? Contact us at: support@ncardia.com www.ncardia.com